



UltraSnap® | Surface ATP Test

UltraSnap® is a user-friendly ATP surface test designed for use with Hygiena® luminometers. This all-in-one test contains a pre-moistened swab bud for optimal sample recovery. A unique liquid-stable reagent provides superior reproducibility and robustness. Hygiena's patented Snap-Valve technology facilitates easy activation.

Benefits:

- Pre-wetted swab maximizes sample recovery and breaks through biofilm
- Unique liquid-stable reagent gives superior sensitivity and reliable results
- 100% recyclable plastic


Sensitivity:

Hygiena System	UltraSnap® Detection Limit				
femtomoles ATP	10 fmol	1 fmol	.1 fmol	.01 fmol	.001 fmol
ATP (mol)	1×10^{-14}	1×10^{-15}	1×10^{-16}	1×10^{-17}	1×10^{-18}
SystemSURE Plus	● ● ● ●				
EnSURE	● ● ● ● ● ●				
EnSURE Touch	● ● ● ● ● ●				

Performance:

UltraSnap's superior performance was validated in a third-party study conducted by a large reference laboratory. Hygiena's system demonstrated superior repeatability, and outstanding performance for linearity, sensitivity, and accuracy over other ATP systems.

Learn more at www.hygiena.com or contact Hygiena for a full copy of the study.



Catalog No.	Description	Quantity
US2020	UltraSnap Surface ATP Test	100



AquaSnap® | Water ATP Tests



AquaSnap® is an easy-to-use ATP test for water samples available in two formats: Total and Free. The test device's unique liquid-stable chemistry, combined with its proprietary test design, offers unrivaled accuracy and reproducibility. The devices are often used for monitoring CIP systems and for industrial water treatment systems such as cooling towers and closed water systems.

Benefits

- Easy to use
- Results in 15 seconds
- Free and Total formats used to reliably and easily indicate microbial contamination in sample
- Proprietary test design enables accurate results with less variation
- Unique liquid-stable reagent provides accuracy and reproducibility
- Reagent is tolerant to temperature abuse and sanitizer effects
- Honey dipper collection tip collects consistent 100 µl sample
- 100% recyclable plastic

Sensitivity

- Extremely sensitive- detects down to 0.1 femtomoles ATP with EnSURE® Touch Monitoring System
- Detects down to 1,000 CFU



Catalog No.	Description	Quantity
AQ-100X	AquaSnap Total Water ATP Test	100
AQ-100FX	AquaSnap Free Water ATP Test	100